

Identifying the Misidentification Problems

Henry Malnick, November 2012



The Pitfalls of Identification

Chadwick, P.R., et al. (1995)
Haemophilus paraphrophilus
infection: a pitfall in laboratory
diagnosis.

J. Infect (1995) 30: 67-69.

Gram-negative catalase negative oxidase positive bacteria

Cardiobacterium hominis

Eikenella corrodens

Kingella denitrificans

Kingella kingae

Suttonella (Kingella) indologenes

Neisseria elongata ss nitroreductans

Aggregatibacter aphrophilus (Haemophilus aphrophilus/paraphrophilus)

Catalase negative, oxidase positive GNR

	ONPG	Lysine	Ornithine	NO ₃ red.	NO ₂ red.	Gas from NO ₂ red.	indole	β-haem
<i>E. corrodens</i>	-	+	+	+	-	-	-	-
<i>N. elongate ss nitroreductans</i>	-	-	-	+	+	-	-	-
<i>Kingella denitrificans</i>	-	-	-	+	+	+	-	-
<i>Kingella kingae</i>	-	-	-	-	-	-	-	+
<i>Cardiobacterium spp. *</i>	-	-	-	-	-	-	+	-
<i>Suttonella indologenes*</i>	-	-	-	-	-	-	+	-
<i>Aggregatibacter spp.</i>	V	V	V	+	-	-	V	V

**Cardiobacterium* produce acid from sorbitol/mannitol, *Suttonella* are negative



Actinobacillus sp.

Aggregatibacter aphrophilus

Haemophilus influenzae

Haemophilus parainfluenzae

Haemophilus sp.


Kingella sp.

Moraxella sp.

Pasteurella multocida,

Pasteurella sp.

Kingella kingae



Misleading clues

- Colony morphology
- Gram stain
- Cell morphology
- Atypical biochemical test results

Gram-positive bacteria submitted as gram-negative bacteria

Arcanobacterium haemolyticum

Aerobic spore bearers (*Brevibacillus*, *Paenibacillus*)

Clostridium tertium

Leifsonia (*Corynebacterium*) *aquaticum*


Exiguobacterium acetylicum

Microbacterium spp.


Propionibacterium acnes


Streptococcus spp.

Capnocytophaga



Leptotrichia spp.








Vancomycin


Gram-positive resistant rods

Lactobacillus, Erysipelothrix


Gram-negative sensitive rods/coccobacilli

Moraxella, Acinetobacter
















**Aerobic spore bearers
identifying as
*Sphingomonas paucimobilis***





A Case of Mistaken Identity











Burkholderia pseudomallei

- **Pseudomonas aeruginosa**
- **Pseudomonas fluorescens**
- **Burkholderia cepacia complex**
- **Chromobacterium violaceum**

Containment level 2 organisms identifying as containment level 3

- **Morganella morganii and Aggregatibacter actinomycetemcomitans identifies as Burkholderia mallei**
- **Chryserobacterium spp. identifies as Francisella tularensis**



Acknowledgement for the Lactobacillus photographs

- **Shila Seaton, Elizabeth Fagan (QAL)**

Acknowledgement for the Bacillus photographs

- **Ijeoma Nwafor, Obioma Mpamugo, Nafeesa Hussain
(Bacillus Reference Unit, HPA Colindale)**

